Jessica Hearn

Winter 2010

# Mein<sub>trup</sub> Fokus

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Hypoxic research helps advance insight into tumour cell metabolism at the Beatson Institute for Cancer Research in Glasgow. Prof Eyal Gottlieb leads a research group studying apoptosis and tumour metabolism. The group is investigating the relationship between cancer, cell metabolism and cell death.

Many laboratory cell culture systems can only expose cells to more than twice the physiological concentration of oxygen found in most cells. The H35 Hypoxystation allows Prof Gottlieb's research to be carried out in a tumour-like hypoxic environment. Many solid tumours consist of large areas of poorly perfused cells leading to areas of low oxygen found throughout the cell mass. As cells are known to behave differently in hypoxic conditions, researchers rely on being able to recreate and maintain a tightly controlled low oxygen environment in order to more accurately reflect the cellular interactions within tumours. The hypoxic tumour regions often produce the most aggressive and therapyresistant cells, causing tumours to spread. Using Don Whitley Scientific's H35 Hypoxystation has allowed researchers to monitor and record consistently low levels of oxygen as well as measuring the alterations in tumour cell metabolism.

Prof Gottlieb studied the role of prolyl hydroxylases (PHD), the enzymes that sense oxygen in cells, in the metabolic adaptation to hypoxia. He observed that by keeping PHD active under hypoxia and hence preventing the hypoxia-mediated metabolic changes tumours would experience a metabolic crisis and undergo cell death. The ability of the Hypoxystation to generate and precisely control an hypoxic environment ensures that the required internal conditions remain constant at all times. Oxygen, carbon dioxide, temperature and humidity levels are continuously displayed on the touchscreen interface. This data, together with chamber pressure, can also be downloaded for further analysis.

Experimenting with the oxygen levels in tumour environment has helped Prof Gottlieb's researchers identify a new approach – direct targeting of PHDs by

addressing hypoxia in general rather than targeting the hypoxia inducible transcription factor (HIF) only. PHDs in most tumours are functionally silenced by the prevailing hypoxia. The research has demonstrated that reactivating these enzymes under hypoxia may be a suitable approach for cancer treatment.

Commenting on the H35 Hypoxystation, Prof Gottlieb commended the design that allows the oxygen content to be changed as and when necessary (in 0.1% increments from 0.1% to 20%). This enables researchers to carry out a comprehensive and more accurate ex vivo analysis of hypoxic tumours that are known to be more aggressive, with high metastatic potential and more resistance to therapy.

Prof Gottlieb's research team continues to investigate the metabolic transformation process that cancer cells undergo in order to survive and grow.

Future studies will continue to focus on the unique metabolic traits of cancer cells and how the survival mechanisms are regulated through the inhibition of PHDs.